


Electricity-Free Incubator

Measurement results



User Information RC-5 Data Logger

Interval 0:15:0 H:M:S

TemperatureMaximum 39.7°C

TemperatureMinimum 36.4°C

Start Time 2018-07-24 15:51:20

End Time 2018-07-26 00:06:20

Temperature Unit °C

Number	Time	Temperature°C
34	2018-07-24 16:06:20	39.7
35	2018-07-24 16:21:20	38.8
36	2018-07-24 16:36:20	38.3
37	2018-07-24 16:51:20	37.8
38	2018-07-24 17:06:20	37.6
39	2018-07-24 17:21:20	37.5
40	2018-07-24 17:36:20	37.3
41	2018-07-24 17:51:20	37.2
42	2018-07-24 18:06:20	37.2
43	2018-07-24 18:21:20	37.2

44	2018-07-24 18:36:20	37.2
45	2018-07-24 18:51:20	37.2
46	2018-07-24 19:06:20	37.3
47	2018-07-24 19:21:20	37.3
48	2018-07-24 19:36:20	37.2
49	2018-07-24 19:51:20	37.3
50	2018-07-24 20:06:20	37.3
51	2018-07-24 20:21:20	37.3
52	2018-07-24 20:36:20	37.4
53	2018-07-24 20:51:20	37.4
54	2018-07-24 21:06:20	37.4
55	2018-07-24 21:21:20	37.4
56	2018-07-24 21:36:20	37.4
57	2018-07-24 21:51:20	37.4
58	2018-07-24 22:06:20	37.4
59	2018-07-24 22:21:20	37.4
60	2018-07-24 22:36:20	37.4
61	2018-07-24 22:51:20	37.4
62	2018-07-24 23:06:20	37.3
63	2018-07-24 23:21:20	37.3
64	2018-07-24 23:36:20	37.3
65	2018-07-24 23:51:20	37.3
66	2018-07-25 00:06:20	37.3
67	2018-07-25 00:21:20	37.3
68	2018-07-25 00:36:20	37.3
69	2018-07-25 00:51:20	37.3
70	2018-07-25 01:06:20	37.3
71	2018-07-25 01:21:20	37.3
72	2018-07-25 01:36:20	37.3
73	2018-07-25 01:51:20	37.3
74	2018-07-25 02:06:20	37.3
75	2018-07-25 02:21:20	37.3
76	2018-07-25 02:36:20	37.3
77	2018-07-25 02:51:20	37.3
78	2018-07-25 03:06:20	37.3
79	2018-07-25 03:21:20	37.3

80	2018-07-25 03:36:20	37.3
81	2018-07-25 03:51:20	37.3
82	2018-07-25 04:06:20	37.3
83	2018-07-25 04:21:20	37.3
84	2018-07-25 04:36:20	37.3
85	2018-07-25 04:51:20	37.3
86	2018-07-25 05:06:20	37.3
87	2018-07-25 05:21:20	37.3
88	2018-07-25 05:36:20	37.3
89	2018-07-25 05:51:20	37.3
90	2018-07-25 06:06:20	37.3
91	2018-07-25 06:21:20	37.3
92	2018-07-25 06:36:20	37.2
93	2018-07-25 06:51:20	37.2
94	2018-07-25 07:06:20	37.2
95	2018-07-25 07:21:20	37.2
96	2018-07-25 07:36:20	37.2
97	2018-07-25 07:51:20	37.2
98	2018-07-25 08:06:20	37.2
99	2018-07-25 08:21:20	37.2
100	2018-07-25 08:36:20	37.2
101	2018-07-25 08:51:20	37.2
102	2018-07-25 09:06:20	37.2
103	2018-07-25 09:21:20	37.2
104	2018-07-25 09:36:20	37.2
105	2018-07-25 09:51:20	37.2
106	2018-07-25 10:06:20	37.2
107	2018-07-25 10:21:20	37.2
108	2018-07-25 10:36:20	37.2
109	2018-07-25 10:51:20	37.2
110	2018-07-25 11:06:20	37.1
111	2018-07-25 11:21:20	37.1
112	2018-07-25 11:36:20	37
113	2018-07-25 11:51:20	37
114	2018-07-25 12:06:20	37
115	2018-07-25 12:21:20	37

116	2018-07-25 12:36:20	37
117	2018-07-25 12:51:20	37
118	2018-07-25 13:06:20	37
119	2018-07-25 13:21:20	37
120	2018-07-25 13:36:20	37
121	2018-07-25 13:51:20	37
122	2018-07-25 14:06:20	37
123	2018-07-25 14:21:20	37
124	2018-07-25 14:36:20	37
125	2018-07-25 14:51:20	36.9
126	2018-07-25 15:06:20	36.9
127	2018-07-25 15:21:20	36.9
128	2018-07-25 15:36:20	36.9
129	2018-07-25 15:51:20	36.9
130	2018-07-25 16:06:20	36.9
131	2018-07-25 16:21:20	36.9
132	2018-07-25 16:36:20	36.9
133	2018-07-25 16:51:20	36.9
134	2018-07-25 17:06:20	36.9
135	2018-07-25 17:21:20	36.9
136	2018-07-25 17:36:20	36.9
137	2018-07-25 17:51:20	36.8
138	2018-07-25 18:06:20	36.8
139	2018-07-25 18:21:20	36.8
140	2018-07-25 18:36:20	36.8
141	2018-07-25 18:51:20	36.8
142	2018-07-25 19:06:20	36.8
143	2018-07-25 19:21:20	36.8
144	2018-07-25 19:36:20	36.8
145	2018-07-25 19:51:20	36.8
146	2018-07-25 20:06:20	36.8
147	2018-07-25 20:21:20	36.8
148	2018-07-25 20:36:20	36.8
149	2018-07-25 20:51:20	36.8
150	2018-07-25 21:06:20	36.7
151	2018-07-25 21:21:20	36.7

152	2018-07-25 21:36:20	36.7
153	2018-07-25 21:51:20	36.7
154	2018-07-25 22:06:20	36.7
155	2018-07-25 22:21:20	36.6
156	2018-07-25 22:36:20	36.7
157	2018-07-25 22:51:20	36.7
158	2018-07-25 23:06:20	36.5
159	2018-07-25 23:21:20	36.5
160	2018-07-25 23:36:20	36.4
161	2018-07-25 23:51:20	36.4
162	2018-07-26 00:06:20	36.4