

**- FOR USE WITH THE MEMBRANE FILTRATION METHOD TO QUANTIFY
BACTERIA IN WATER SAMPLES -**



GENERAL INFORMATION

Cellulose Nitrate membranes for microbiological water analysis by membrane filtration according to ISO 9308-1: 2014.

ADVANTAGES

- Ready-to-use
- Pre-sterilized
- Grid print for easy counting
- Easy to transport and store



SPECIFICATIONS

- Material: Cellulose Nitrate
- Diameter: 47 mm
- Pore width: 0,45 µm
- Shelf life: 4 years after production
- Grid print
- Individually and sterile packed